



Received on 16 September, 2015; received in revised form, 05 December, 2015; accepted, 05 January, 2016; published 01 March, 2016

DEVELOPMENT AND VALIDATION OF ANALYTICAL METHOD FOR DETERMINATION OF MICAFUNGIN AND ITS RELATED SUBSTANCES IN BULK BY RP-UPLC

Smita Joshi ^{*1}, Falguni Majmudar ² and Nirav Vyas ³

Department of Pharmaceutical Chemistry ¹, K. B. Raval College of Pharmacy, Gandhinagar, Gujarat, India
Department of Pharmacology ², Smt. N.H.L. Municipal Medical College, Ahmedabad, Gujarat, India
Analytical Research development ³, Amneal Pharmaceuticals Pvt. Ltd., Ahmedabad, Gujarat, India

Key words:

Micafungin sodium, RP-UPLC, related substances, impurities

Correspondence to Author:

Smita Joshi

Department of Pharmaceutical Chemistry, K. B. Raval College of Pharmacy, Gandhinagar, Gujarat, India.

Email: smita_talaviya85@yahoo.com

ABSTRACT: A gradient reverse phase Ultra Performance Liquid Chromatographic (RP-UPLC) method was developed and validated for determination of Micafungin sodium and its synthetic impurities. The successful separation of Micafungin sodium and its synthetic impurities was achieved using Phenomenax Aeris peptide XB C18 (150×2.1 mm i.d., 1.7 μ particle size) column maintained at 45 °C temperature with mobile phase consisting 0.01 M phosphate buffer pH 2.9 and acetonitrile in a gradient programme. The mobile phase flow rate was 0.3 ml/min and the detection wavelength was 279 nm. The retention time in developed method is comparatively less than the reported HPLC methods for determination of Micafungin sodium offering less time consuming and fast analytical method. The developed RP-UPLC method was validated according to ICH guidelines with respect to linearity, accuracy, precision, specificity and robustness and also the LOD and LOQ values were determined.

INTRODUCTION: Micafungin is a semi synthetic compound, belonging to a newer class of antifungal agents, echinocandins lipopeptides. It is synthesized by chemical modification of a fermentation product from *Coleophoma empetri* ¹. Micafungin selectively inhibit 1,3-β-D-glucan, which is required for fungal cell wall synthesis. It has been approved for the treatment of esophageal candidiasis and for prophylaxis of *Candida* and *Aspergillus* infections in patients undergoing hematopoietic stem cell transplantation ^{2,3}.

The chemical structure of Micafungin sodium, impurity A, impurity B, impurity C and impurity D are shown in figure 1, the molecular weight of the compound is 1292.26 gm/mol and the empirical formula is C₅₆H₇₀N₉O₂₃S ⁴. The drug was first launched in Japan in December 2002 and then also approved by US food and drug administration in March 2005 ^{5,6}.


The quality and safety of drug product can be significantly affected by the presence of impurities and therefore it is important to study the impurity profile of drug substance ⁷. Ultra performance liquid chromatography (UPLC) is relatively novel technique with new possibilities in liquid chromatography, focusing on decrease of analysis time and solvent consumption ⁸. UPLC systems can withstand high system back pressure and comprising special analytical columns with sub 2

QUICK RESPONSE CODE	DOI: 10.13040/IJPSR.0975-8232.7(3).1211-18
	Article can be accessed online on: www.ijpsr.com
DOI link: http://dx.doi.org/10.13040/IJPSR.0975-8232.7(3).1211-18	


μm particles. UPLC system allows shortening analysis time up to nine times compared to conventional HPLC system⁹.

Few HPLC methods have been reported in the literature for quantification of Micafungin in plasma¹⁰⁻¹³. Shengsheng et al reported stability indicating HPLC method for determination of


Micafungin and related substances¹⁴. To the best of our knowledge no UPLC method has been reported for the estimation of Micafungin with its related substances. The aim of this work was to develop a rapid, simple, precise, accurate and validated UPLC method for determination of Micafungin and its related substances.




1(a)




1(b)



1(c)



1(d)



1(e)

FIG.1: MICA FUNGIN SODIUM (1a), IMPURITY A - MICA FUNGIN RING OPENING IMPURITY (1b), IMPURITY B - MICA FUNGIN C-26 DESMETHYL IMPURITY (1c), IMPURITY C - MICA FUNGIN C-21 S-IMPURITY (1d), IMPURITY D - MICA FUNGIN DESULFATE IMPURITY (1e).

Experimental:**Chemicals and Instrumentation:**

Samples of Micafungin and its four impurities were received as gift sample from Amneal Pharmaceuticals. Specification levels for Impurity A, Impurity B, Impurity C and Impurity D were $\leq 0.5\%$, $\leq 0.15\%$, $\leq 0.15\%$ and $\leq 0.15\%$ respectively. The UPLC systems, (Shimadzu, Nexera system) consist of autoinjector with PDA detector. The analytical column, Phenomenex Aeris peptide XB C18 (150*2.1 mm i.d., 1.7 μ particle size), was operated at 45 °C temperature. The gradient Programme consists of 0.01 M phosphate buffer, pH 2.9 was adjusted with 85 % orthophosphoric acid solution in water and acetonitrile (HPLC grade).

Stock and working standard solutions:

Accurately weighed quantity of Micafungin, 50 mg was transferred into separate 50 ml volumetric flask, dissolved and diluted up to mark with diluent (Phosphate buffer: Acetonitrile 95:5 v/v) to prepare 1000 μ g/ml of stock solution. A stock solution of impurities (Mixture of imp-A, imp-B, imp-C and imp-D) were also prepared in the diluent.

Method Validation:

After method development, validation of the chromatographic method for Micafungin was performed in accordance with ICH guidelines¹⁵.

Linearity and Range:

Linearity is the ability of a method to elicit test results that are directly proportional to analyte concentration within a given range. Range is the interval between the upper and lower levels of analyte that have been demonstrated to be determined with precision, accuracy and linearity. Linearity of the method was checked at seven different concentrations. Linearity solutions of 50, 100, 150, 200, 250, 300 and 350 μ g/ml were prepared by transferring 0.5, 1, 1.5, 2, 2.5, 3 and 3.5 ml of stock solution to 10 ml of volumetric flask and volume was adjusted with diluent. Linearity solutions for impurity were prepared at 150, 100, 75, 50 and 25 % specification level for each impurity.

Specificity:

Specificity is the ability of the method to measure the analyte response in the presence of its potential

impurities. The specificity of the developed method for Micafungin was carried out in the presence of its impurities namely, imp-A, imp-B, imp-C and imp-D.

Accuracy:

Accuracy of analytical method is the closeness of test results to the true value. Accuracy was determined over the range from 50 % as lowest sample concentration to 150 % as highest sample concentration. Triplicate preparations for each level were prepared and injected in presence of spiked impurities to their 100 % specification level. The mean recovery should be within 100 ± 2 % at each concentration over the range of 50-150 % of concentration. Accuracy was demonstrated at 50 %, 100 % and 150 % level with respect to specification limit spiked with API for each impurity. Recovery was determined in triplicates and reported.

Precision:

The precision of the method was evaluated by carrying out six independent assays of Micafungin test samples and calculating the % relative standard deviation (RSD) of the assay. Precision with related substances was determined by spiking the impurities with Micafungin sodium (300 μ g/ml) at 100 % specification level of each impurity (six replicates). % RSD of area of each imp-A, imp-B, imp-C and imp-D was calculated. Intraday and interday precision was carried out and % RSD was determined.

Limit of Detection and Limit of Quantification:

The limit of detection and limit of quantification were evaluated by serial dilutions of Micafungin sodium, Impurity A, Impurity B, Impurity C and Impurity D stock solutions in order to obtain signal to noise ratio of 3:1 for LOD and 10:1 for LOQ.

Robustness study:

Robustness study was carrying out by changing the minor parameters of the chromatographic conditions. The assay was performed by change in the flow rate, column temperature and minor change in pH by taking three replicates and % RSD was calculate. Resolution was determined between the main peak and closest impurity i.e. impurity B.

RESULT AND DISCUSSION:**Optimization of Chromatographic conditions:**

The main objective of the chromatographic method development was to separate Micafungin sodium from the imp-A, imp-B, imp-C and imp-D. Several trials were carried out for accurate and precise method development and impurities were coeluted. After using several columns and buffers, suitable column chemistry and good peak shape were

obtained with Phenomenax Aeris peptide XB C18 (150×2.1 mm i.d., 1.7 μ particle size), column temperature was adjusted at 45 °C, with gradient mobile phase system consisting Mobile phase A as 0.01M phosphate buffer and Mobile phase B as Acetonitrile. The best separation of Micafungin from the impurities was achieved with a gradient program as follows:



TABLE 1: GRADIENT PROGRAM IN OPTIMIZED CHROMATOGRAPHIC CONDITIONS

Time (minutes)	Flow rate (ml/min)	% of Mobile phase A	% of Mobile phase B
0	0.3	80	20
5	0.3	80	20
10	0.3	60	40
15	0.3	20	80
20	0.3	80	20
25	0.3	80	20

Where, Mobile phase A is 0.01M phosphate buffer and Mobile phase B is Acetonitrile

The flow rate was 0.3 ml/min. The injection volume was 5 μ l. Depending on peak absorption maxima and peak purity index, 279 nm was selected for detection. The sample temperature was kept at 5 °C. In optimized chromatographic conditions, Micafungin, imp-A, imp-B, imp-C and imp-D were separated with good resolution, typical

retention time were 15.14, 14.41, 14.92, 16.10 and 16.89 min respectively (**Fig. 2**). The system suitability results are given in Table and the developed method was found to be specific for Micafungin and its known impurities namely imp-A, imp-B, imp-C and imp-D.

**FIG.2: CHROMATOGRAM OF MICA FUNGIN SODIUM (RT-15.14 Min), IMPURITY A (Rt-14.41), IMPURITY B (Rt-14.92), IMPURITY C (RT-16.10) AND IMPURITY D (Rt-16.89) IN THE OPTIMIZED CHROMATOGRAPHIC CONDITIONS****FIG. 3: CHROMATOGRAM OF DILUENT (BLANK) IN THE OPTIMIZED CHROMATOGRAPHIC CONDITIONS**

Linearity and Range: Linearity of the method was evaluated at seven different concentrations, at 50, 100, 150, 200, 250, 300 and 350 µg/ml. The peak area was recorded and plotted versus standard concentrations. Results shown that the method was linear as straight line equation, $y = 24242.71x - 11003.7$, was obtained, where y indicates peak area and x indicates concentration of Micafungin over the specified range with R² of 0.9999 (Table 2 and Fig. 4). For imp-A, imp-B, imp-C and imp-D, 0.9995, 0.9995, 0.9999 and 0.9996 R² value was obtained.

TABLE 2: LINEARITY DATA FOR MICA FUNGIN SODIUM

Concentration (µg/ml)	Mean area	%RSD
50	1208143	0.179
100	2407086	0.258
150	3633700	0.187
200	4829775	0.156
250	6040062	0.059
300	7259917	0.075
350	8484064	0.092




FIG. 4: CALIBRATION CURVE OF MICA FUNGIN SODIUM (50-350 µg/ml)

Specificity

The developed analytical method was found to be specific as the peak purity index for Micafungin sodium was 999.11 in presence of all the spiked impurities i.e. Impurity A, Impurity B, Impurity C and Impurity D; and also the resolution was more than 1.5 between each peaks.

Accuracy

% Recovery for Micafungin sodium at 50 %, 100 % and 150 % concentration level was found to be 100.03, 100.08 and 99.83 respectively. % Recovery for Impurity A, Impurity B, Impurity C and Impurity D at 50 %, 100 % and 150 % specification level was found to be between 98-102 %, within the acceptance criteria, so the method is accurate for determination of Micafungin sodium and its impurities (Table 3)

TABLE 3: ACCURACY RESULTS FOR MICA FUNGIN SODIUM AND ITS IMPURITIES

Name of substances	Concentration level / specification level	% Mean Recovery ± S.D
Micafungin sodium	50 %	100.03±0.39
	100 %	100.08±0.17
	150 %	99.83±0.17
Impurity A	50 %	100.19±0.95
	100 %	99.91±0.47
	150 %	99.95±0.33
Impurity B	50 %	100.27±1.04
	100 %	99.63±0.78
	150 %	100.13±0.37
Impurity C	50 %	98.95±0.51
	100 %	98.33±0.35
	150 %	99.82±1.45
Impurity D	50 %	98.56±0.69
	100 %	98.88±0.88
	150 %	99.05±1.4

Precision:

The relative standard deviation was found to be 0.083, 0.575, 0.817, 0.954 and 0.809 for Micafungin sodium, imp-A, imp-B, imp-C and imp-D respectively in the repeatability study. % RSD of peak area is less than 2.0 %, for all impurities and Micafungin sodium in the repeatability study (Table 4). Intraday and interday precision also exhibited % RSD less than 2 (Table 5). So the developed method is precise for its use.

TABLE 4: PRECISION (REPEATABILITY) DATA FOR MICA FUNGIN SODIUM, IMPURITY A, IMPURITY B, IMPURITY C AND IMPURITY D.

Result	Micafungin sodium	Impurity A	Impurity B	Impurity C	Impurity D
Mean Area	7259209	19587.170	7022.667	1975.833	2378.333
% RSD	0.083	0.575	0.817	0.954	0.809

TABLE 5: INTRADAY AND INTERDAY PRECISION FOR MICAFUNGIN SODIUM

Intraday precision			Interday precision		
Conc. (µg/ml)	Area Mean	% R.S.D	Conc. (µg/ml)	Area Mean	% R.S.D
200	4835775	0.138	200	4831040	0.169
250	6040462	0.086	250	6042476	0.119
300	7259727	0.119	300	7263306	0.112

Limit of Detection and Limit Quantitative

The LOD and LOQ values for Micafungin sodium was found to be 0.02 and 0.05 µg/ml. LOQ value was precised by six replicate injections with % RSD of 0.364 and checked for linear response with respect to other linearity levels by extending linearity curve upto LOQ level, with R² of 0.9999. The chromatogram of Micafungin sodium at LOQ level is shown in Fig. 5.




FIG.5: CHROMATOGRAM OF MICAFUNGIN SODIUM AT LOQ (0.05 µg/ml) LEVEL

Robustness study:

With the minor modifications in the chromatographic conditions i.e. Flow rate, pH of mobile phase and temperature of column, the % assay was not affected and % RSD was found to be less than 2 % and the resolution between Micafungin sodium and Impurity B was not significantly affected (Table 6). Also with these minor modifications, System suitability parameters were also not affected and hence, the developed analytical method is robust.

TABLE 6: ROBUSTNESS STUDY

Parameter	Assay results for Micafungin sodium		Resolution between Micafungin sodium and Impurity B
	% Assay±SD (n=3)	% RSD	
Standard	99.63±0.20	0.204	1.51
Flow rate (ml/min)	100.01±0.29	0.298	1.48
pH of mobile phase (% v/v)	99.65±0.78	0.783	1.50
Column	99.88±0.46	0.469	1.52
temperature (°C)	99.68±0.31	0.311	1.50
	100.02±0.88	0.883	1.53
	99.24±0.08	0.086	1.51

System suitability:

System suitability parameters were calculated at the start of study of each validation parameter. The values of system suitability results obtained during the entire study are recorded. System Suitability values are from the first injection of six replicates of standard and are derived from Empower 3 software (Table 7).

TABLE 7: SYSTEM SUITABILITY DATA

Compound (n=6)	Retention time (minutes)	USP tailing	USP resolution	USP plate count
Impurity A	14.411	1.13	-	106883
Impurity B	14.920	0.93	2.63	77306
Micafungin sodium	15.143	1.41	1.52	146955
Impurity C	16.100	1.67	6.07	159026
Impurity D	16.895	0.89	4.34	131548

CONCLUSION: A novel, reverse phase liquid chromatographic method has been developed and validated for the estimation of Micafungin sodium and its impurities namely impurity A, impurity B, impurity C and impurity D with a very recent and advanced UPLC method. The proposed method is found to be simple, accurate, precise, sensitive, specific and robust. Hence, it can be successfully used for the routine analysis of Micafungin sodium in pharmaceutical dosage forms and can be also successfully utilized for estimation of Micafungin sodium and its related substances in bulk and in dosage form.

ACKNOWLEDGEMENT: The authors are thankful to Amneal Pharmaceuticals Pvt. Ltd., for providing the use of instruments and materials.

REFERENCES:

1. Tawara S, Ikeda F, Maki K, Morishita Y, Otomo K, Teratani N, Goto T, Tomishima M, Ohki H, Yamada A: In vitro activities of a new lipopeptide antifungal agent, FK463, against a variety of clinically important fungi. Antimicrob. Agentschemother. 2000; 44:57-62.
2. Scott LJ: Micafungin: a review of its use in the prophylaxis and treatment of invasive Candida infection. Drugs 2012;72(16):2141-65

3. Hatano K, Morishita Y, Nakai T, Ikeda F: Antifungal mechanism of FK463 against *Candida albicans* and *Aspergillus fumigatus*. *J. Antibiot.* 2002;55:219-222
4. Veherschild JJ, Cornely OA: Micafungin sodium, the second of the echinocandin class of antifungals: Theory and practice. *Future Microbiol.* 2006;1:161-170
5. Carver PL: Micafungin. *Ann. Pharmacother.* 2004;38:1707-21
6. Zaas AK, Steinbach WJ: Micafungin: The US perspective. *Expert Rev. Anti-infect. Ther.* 2005;3:183-190
7. Bakshi M, Singh S: Development of validated stability-indicating assay methods-Critical review. *J. Pharm. Biomed. Anal.* 2002;28:1011-1040
8. Novakova L, Matysova L, Solich P: Advantages of application of UPLC in pharmaceutical analysis. *Talanta* 2006;68:908-918
9. Michael ES: UPLC: An introduction and review. *J Liq. Chromatogr. Rel. techno.* 2005;28:1253-1263
10. Noda K, Yokoyama T, Ankiyama Y, Tadano K: A liquid chromatographic method for quantifying unbound micafungin in human plasma and the relationship between the concentrations of unbound and total micafungin in plasma. *Yakugaku Zasshi* 2013;133:397-404
11. Nakagawa S., Kuwabara N., Kobayashi H., Shimoeda S., Ohta S., Yamato S: Simple column-switching HPLC method for determining levels of the antifungal agent micafungin in human plasma and application to patient samples. *Biomed. Chromatogr.* 2013;27:551-555
12. Martens-Lobenhoffer J., Rupprecht V., Bode-Boger S.M: Determination of micafungin and anidulafungin in human plasma: UV- or mass spectrometric quantification? *J. Chromatogr. B.* 2011;879:2051-2056
13. Uranishi H., Nakamura M., Nakamura H., Ikeda Y., Otsuka M., Kato Z., Tsuchiya T: Direct-injection HPLC method of measuring micafungin in human plasma using a novel hydrophobic/hydrophilic hybrid ODS column. *J. Chromatogr. B.* 2011;879:1029-1032
14. Shengsheng Z, Xiang M, Xin S, Yongwei L, Zuyue S: Development and Validation of a Stability-Indicating High Performance Liquid Chromatographic (HPLC) method for the determination of related Substances of Micafungin Sodium in drug substances. *Int J Mol Sci* 2013; 14(11):21202-21214
15. ICH, Q2 (R1), Validation of Analytical Procedures: Text and Methodology: 2005

How to cite this article:

Joshi S, Majmudar F and Vyas N: Development and Validation of Analytical Method for Determination of Micafungin and Its Related Substances in Bulk by RP-UPLC. *Int J Pharm Sci Res* 2016; 7(3): 1211-18. doi: 10.13040/IJPSR.0975-8232.7(3).1211-18.

All © 2013 are reserved by International Journal of Pharmaceutical Sciences and Research. This Journal licensed under a Creative Commons Attribution-NonCommercial-ShareAlike 3.0 Unported License.

This article can be downloaded to **ANDROID OS** based mobile. Scan QR Code using Code/Bar Scanner from your mobile. (Scanners are available on Google Playstore)