



Received on 13 May, 2017; received in revised form, 03 January, 2018; accepted, 28 January, 2018; published 01 February, 2018

DEVELOPMENT AND VALIDATION OF AN RP-HPLC METHOD FOR QUANTITATIVE ESTIMATION OF ILOPERIDONE IN RAT PLASMA

V. Dubey and T. R. Saini*

Industrial Pharmacy Research Lab, Department of Pharmacy, Shri G. S. Institute of Technology and Science, Indore - 452003, Madhya Pradesh, India.

Keywords:

Iloperidone, RP-HPLC,
Method validation, PDA detector,
Protein precipitation, Rat plasma

Correspondence to Author:

Dr. T. R. Saini

Industrial Pharmacy Research Lab,
Department of Pharmacy, Shri G. S.
Institute of Technology and Science,
Indore - 452003, Madhya Pradesh,
India.


E-mail: profsaini@gmail.com

ABSTRACT: A simple and fast isocratic reverse phase- high performance liquid chromatography (RP-HPLC) method was developed for estimation of iloperidone from rat plasma employing a photo diode array (PDA) detector. The plasma samples were prepared by protein precipitation using acetonitrile as extracting solvent. The chromatographic separation was performed on a BDS Hypersil Phenyl column (250 mm × 4.6 mm; 5 μm) using a mixture of 0.01 M sodium dihydrogen phosphate buffer (pH 6.0) and acetonitrile (60:40, v/v) as mobile phase with a flow rate of 1.5 ml/min, retention time 7.04 min, and detection wavelength 276 nm. The method was validated for linearity, sensitivity, selectivity, accuracy, precision, and stability according to USFDA guidelines. The calibration curve exhibited linearity in the concentration range of 50–1000 ng/ml ($r^2=0.998$) and the lower limit of quantification (LLOQ) and limit of detection (LOD) of the method were 50 ng/ml and 30 ng/ml respectively. The chromatographic runs were selective without any interfering peaks at the retention time of iloperidone. The accuracy and precision of the method were within USFDA specified range. The plasma samples showed stability over a period of 24 h at room temperature, for two months at - 20 °C, and three freeze/thaw cycles.

INTRODUCTION: Iloperidone is an atypical antipsychotic drug approved in the year 2009 by USFDA for the treatment of schizophrenia^{1,2}. It is quite effective in dealing with both positive and negative symptoms of schizophrenia and has many advantages over other atypical antipsychotic agents. Its tolerability in patients is exceedingly good as it induces very less extrapyramidal symptoms³. Currently it is marketed as tablet dosage form in the strengths of 1 - 12 mg.

Several HPLC methods have been reported for quantification of iloperidone in bulk drug and tablet dosage forms^{4,5,6,7}. For its estimation in human plasma, LCMS methods have been developed^{8,9,10,11,12}. The LCMS based methods however, are not only expensive but the non-availability of this instrumental facility in the conventional research labs remains a constraint as compared to HPLC methods which are most widely employed in the drug analysis.

Although, an HPLC method has been reported for simultaneous determination of iloperidone in presence of idebenone in rat plasma using UV detector with 8.108 min retention time¹³, but the HPLC method employing a PDA detector has not been developed for the estimation of iloperidone in

<p>QUICK RESPONSE CODE</p> 	<p>DOI: 10.13040/IJPSR.0975-8232.9(2).790-94</p> <hr/> <p>Article can be accessed online on: www.ijpsr.com</p>
<p>DOI link: http://dx.doi.org/10.13040/IJPSR.0975-8232.9(2).790-94</p>	

biological samples. The PDA detector facilitates faster analyte quantification at multiple wavelengths simultaneously and also possesses a higher wavelength precision, sensitivity, reliability, and ruggedness in comparison to UV detector^{14,15}.

In the present study a simple, fast, sensitive, and reproducible HPLC method was developed for quantitative estimation of iloperidone in rat plasma using PDA detector.

MATERIALS AND METHODS: Iloperidone was supplied as gift sample by M/S Sun pharmaceutical industries, Mumbai, India. HPLC grade acetonitrile was purchased from Merck Specialities, Mumbai. Ultrapure water was obtained from Millipore purification system (Direct-Q[®]3-R, Merck Millipore, India) of our lab. The membrane filter (0.45 µm) was purchased from Millipore, India. All other chemicals were of analytical grade and purchased from the market.

Sample Preparation: The blood sample was collected from lateral tail vein of albino rats in centrifuge tubes containing disodium edetate solution as anticoagulant and centrifuged at 5000 rpm for 10 min at 4 °C (Eppendorf 5415R Centrifuge, Eppendorf, Germany) to separate the plasma. Iloperidone stock solution of 1 mg/ml concentration was prepared by dissolving 100 mg of iloperidone in 100 ml of acetonitrile. Standard solutions of different concentrations of iloperidone were prepared from this stock solution by diluting with acetonitrile. Standard solution of iloperidone in rat plasma was prepared by spiking 100 µl of the rat plasma with 100 µl of iloperidone standard solution, adding 200 µl of acetonitrile and then mixing on a vortex mixer (CM-101 Plus, Remi Lab World, India) for 10 min to completely precipitate plasma proteins. This mixture was centrifuged at 5000 rpm for 10 min at 4 °C (Eppendorf 5415R Centrifuge, Eppendorf, Germany) and the supernatant was collected in another tube and subjected to chromatographic analysis.

Chromatographic Conditions: The chromatographic analysis was performed on a Waters reverse phase HPLC system (Waters Corporation, USA) equipped with a Waters 2996 photodiode array detector, Waters 600 pump, and Empower software using a BDS Hypersil Phenyl column (4.6

mm × 250 mm; particle size 5 µm; Thermo Fisher Scientific Inc., USA) The mobile phase was consisted of a mixture of 0.01 M sodium dihydrogen phosphate buffer (pH 6.0) and acetonitrile (60:40, v/v) with a flow rate of 1.5 ml/min. The injection volume was 20 µl using a fixed loop, total run time 8.5 min, and the detector wavelength 276 nm.

Method Validation: The validation of the developed analytical method was done in accordance with the USFDA guidelines for bioanalytical method validation¹⁶.

Linearity: The standard solutions of iloperidone in rat plasma were prepared in the concentration range of 50 - 1000 ng/ml (n=6). Each sample was subjected to chromatographic analysis starting from lower to higher concentration. A standard curve was constructed by plotting chromatographic peak area against the drug concentration and its linearity was statistically confirmed.

Sensitivity: The minimum concentration of analyte which can be accurately and precisely measured was calculated using signal to noise ratio. Serially diluted standard solutions of iloperidone were spiked into the rat plasma and analyzed (n=6). The concentrations yielding three times and five times more intense chromatographic response as compared to the baseline noise were measured as the limit of detection (LOD) and lower limit of quantification (LLOQ) of iloperidone respectively.

Selectivity: The blank plasma samples obtained from six different rats were separately analyzed (n=6) and the resulting chromatograms were compared with those obtained from the analysis of standard solutions of iloperidone in rat plasma at the lower limit of quantification in order to find out any possible interference of the endogenous substances present in the plasma at the retention time of the peak.

Accuracy: Recovery method was used to determine accuracy of the developed RP-HPLC method. For this three concentrations of standard solution of iloperidone in rat plasma, *i.e.*, 100 ng/ml, 400 ng/ml, and 1000 ng/ml (nominal drug concentrations) respectively were estimated (n=6) and the mean recovery was calculated. The

accuracy was reported as percent relative error (%RE) of the recovered amount.

Precision: Precision of the developed method was assessed by analyzing the standard solution of iloperidone in rat plasma (n=6) in the strengths of 100 ng/ml, 400 ng/ml, and 1000 ng/ml (nominal drug concentrations) and measuring peak areas of their chromatograms. The within-run precision, *i.e.*, repeatability of the method within a batch was determined by estimating the nominal drug concentrations in the same run on the same day, while the between-run precision, *i.e.*, repeatability of the method over time was determined by estimating them on three different days. The precision was measured as the percent coefficient of variation.

Stability: The stability of standard solutions of iloperidone in rat plasma was determined (n=3) both at a lower concentration (100 ng/ml) and a higher concentration (1000 ng/ml) at room temperature for twenty four hours (bench top condition), at - 20 °C for two months (long-term storage), and after subjecting to three freeze/thaw cycles alternating - 20 °C and room temperature to mimic the sample handling conditions.

RESULTS AND DISCUSSION:

Optimization of Chromatographic Conditions: Different types of columns and mobile phases were evaluated for obtaining a selective and efficient separation with a symmetric peak and optimum retention time. Three types of columns namely, Waters C - 8 Column, Thermo BDS Hypersil C-18 column, and Thermo BDS Hypersil Phenyl column were studied. Since, iloperidone is practically insoluble in water and freely soluble in acetonitrile¹⁷, the mixtures of phosphate buffer of different pH and acetonitrile in different ratios were employed as mobile phases.

The BDS Hypersil Phenyl column (4.6 mm × 250 mm, particle size 5 µm) provided the best selectivity and sharp peaks with maximum signal to noise (S/N) ratio and minimum retention time as compared to Waters C-8 and Thermo BDS Hypersil C-18 columns which resulted in significant peak tailing and higher retention times (> 12 min). The BDS (base deactivated silica) columns are more selective in nature because they

are end-capped for reduced tailing of the peak and attached with phenyl group to provide better selectivity. The pH of the buffer used in mobile phase was found to affect the peak shape such that the extremes of pH towards both acidic and basic sides resulted in significant peak tailing. At pH 6.0 the sharpest peak with least tailing was observed.

The mixture of 0.01 M sodium dihydrogen phosphate buffer (pH 6.0) and acetonitrile in 60:40 ratio was selected as the mobile phase because it produced a sharp and resolved peak shape with only 40% of the organic solvent. The spectral peak of iloperidone in the non-interfering zone was most prominent at 276 nm and therefore it was selected as the detection wavelength. The retention time for elution of iloperidone was 7.04 min with a total run time of 8.5 min (**Fig. 1**).

Extraction Procedure: The protein precipitation method using acetonitrile as solvent was employed for extraction of iloperidone from plasma because it is a simple method and involves lesser number of steps¹⁸. Acetonitrile was selected as solvent because iloperidone is freely soluble in it and it also has a good plasma protein precipitating property.

Method Validation:

Linearity and Calibration Curve: The calibration curve constructed between the chromatographic peak area of iloperidone against its concentration was linear over the concentration range of 50–1000 ng/ml with a correlation coefficient (r^2) of 0.998. This confirmed the linearity of the developed chromatographic method.

Sensitivity: The least concentration of the analyte (LLOQ) measured with the acceptable accuracy and precision from the standard solution of iloperidone in rat plasma was found to be 50 ng/ml and the limit of detection (LOD), *i.e.*, the least quantity of analyte which gives a peak area of three times that of baseline noise (signal-to noise ratio of 3:1) was detected as 30 ng/ml.

Selectivity: The chromatogram of spiked samples showed sharp peak of iloperidone at its retention time. The chromatogram of the blank plasma showed only few peaks between 1 and 4 minutes of detection time and no significant signals were seen afterwards (**Fig. 1**).

The absence of any peak at the retention time of analyte showed that there is no significant interference of rat plasma at the retention time of iloperidone in spiked samples. This established the selectivity of the developed method for analysis of iloperidone.




FIG. 1: HPLC CHROMATOGRAMS OF BLANK PLASMA [A] AND PLASMA SPIKED WITH ILOPERIDONE [B]

Accuracy: As shown in **Table 1**, the intra-day and inter-day percent relative error (%RE) of the recovered amount of iloperidone from its standard solutions in rat plasma was found to be 4.64%, -5.67% and 3.68% and 5.80%, 2.96% and 4.20% respectively for 100 ng/ml, 400 ng/ml, and 1000 ng/ml concentrations. Since the recovery of drug was in the range of $\pm 6\%$ of the nominal drug concentrations, *i.e.*, well within the limits of $\pm 15\%$ as prescribed by USFDA¹⁶, it confirmed that the developed method can accurately determine iloperidone in rat plasma.

Precision: Precision was determined by analyzing three drug concentrations of standard solutions of iloperidone in rat plasma, *i.e.*, 100ng/ml, 400ng/ml, and 1000ng/ml. As shown in **Table 1**, the % coefficient of variation (%CV) for within-run or intra-day precision was between 2.02% and 8.77% and that for between-run or inter-day precision on three different days was between 2.84% and 6.31%. Since the values of both intra-day and inter-day precision were within the accepted limits of 15% CV as prescribed by USFDA,¹⁶ it confirmed that the developed method is a precise method for analysis of iloperidone in rat plasma.

TABLE 1: ACCURACY AND PRECISION DATA

Nominal concentration (ng/ml)	Mean measured concentration (ng/ml)	Accuracy RE (%)	Precision CV (%)
Within-run/intra-day			
100	104.64	4.64	4.65
400	377.33	-5.67	2.02
1000	1036.81	3.68	8.77
Between-run/inter-day			
100	105.80	5.80	3.39
400	411.84	2.96	2.84
1000	1042.02	4.20	6.31

Stability: The results presented in **Table 2** show that the percent nominal concentration of standard solutions of iloperidone in rat plasma subjected to stability study was found to be within the range of 4% of the nominal drug concentrations under all

the stability conditions which was much below the 15% limit prescribed by USFDA¹⁶. This confirmed that the standard solutions of iloperidone in rat plasma were stable.

TABLE 2: RESULTS OF STABILITY STUDY

Stability study type	Stability condition	Nominal concentration (ng/ml)	Concentration after study* (ng/ml)	RE (%)	% Nominal concentration	CV (%)
Bench-top	24 hours at RT	100	102.83	2.83	102.83	2.91
		1000	996.54	-0.35	99.65	2.17
Long-term	2 months at -20 °C	100	98.31	-1.69	98.31	3.28
		1000	970.84	-2.92	97.08	7.84
Freeze-thaw	3 cycles	100	103.86	3.86	103.86	4.41
		1000	982.87	-1.71	98.29	5.12

*Mean (n = 3)

CONCLUSION: A simple, fast, sensitive, accurate, and precise RP-HPLC method using PDA detector was developed for estimation of iloperidone in rat plasma with a retention time of 7.04 min and total run time of 8.5 min. The mobile phase used in the method consisted of only 40% organic solvent. The method was linear over a concentration range of 50 - 1000 ng/ml ($r^2 = 0.998$) with an LLOQ of 50 ng/ml and demonstrated good accuracy and precision. The iloperidone was stable in rat plasma samples for 24 h at room temperature, two months at - 20 °C, and under three freeze/thaw cycles.

ACKNOWLEDGEMENT: Authors are thankful to M/S Sun Pharmaceutical Industries Ltd., Mumbai, India for providing the gift sample of iloperidone.

CONFLICT OF INTEREST: There are no conflicts of interest.

REFERENCES:

- Weiden PJ: Iloperidone for the treatment of schizophrenia: An updated clinical review. *Clin Schizophr Relat Psychoses* 2012; 6: 34 - 44.
- Citrome L: Iloperidone: Chemistry, pharmacodynamics, pharmacokinetics and metabolism, clinical efficacy, safety and tolerability, regulatory affairs, and an opinion. *Expert Opin Drug Metab Toxicol* 2010; 6: 1551-64.
- Cutler AJ, Kalali AH, Weiden PJ, Hamilton J and Wolfgang CD: Four-week, double-blind, placebo- and ziprasidone-controlled trial of iloperidone in patients with acute exacerbations of schizophrenia. *J Clin Psychopharmacol* 2008; 28: S20-8.
- Mandpe LP and Pokharkar VB: Stress degradation studies on Iloperidone and development of a stability indicating HPLC method for bulk drug and pharmaceutical dosage form. *Der Chemica Sinica* 2011; 2: 230-9.
- Chhalotiya UK, Bhatt KK, Shah DA and Patel JR: Liquid Chromatographic Method for the Quantification of Antipsychotic Agent Iloperidone in Pharmaceutical Formulation. *ISRN Analytical Chemistry* 2012; 2012: 6.
- Landge SB, Jadhav SA, Vishwambar SP, Solanki PV, Bembalkar SR and Mathad VT: Development and Validation of RP-UPLC Method for the Determination of Iloperidone, Its Related Compounds and Degradation Products in Bulk and Dosage Form. *Am J Anal Chem* 2014; 5: 13.
- Atila Karaca S and Uğur DY: Determination of Iloperidone in Pharmaceuticals by Validated High-Performance Liquid Chromatography. *Analytical Letters* 2017: 1-9.
- Mutlib AE and Strupczewski JT: Picogram determination of iloperidone in human plasma by solid-phase extraction and by high-performance liquid chromatography-selected-ion monitoring electrospray mass spectrometry. *J Chromatogr B Biomed Appl* 1995; 669: 237-46.
- Parekh JM, Sanyal M, Yadav M and Shrivastav PS: Stable-isotope dilution LC-MS/MS assay for determination of iloperidone and its two major metabolites, P 88 and P 95, in human plasma: application to a bioequivalence study. *Bioanalysis* 2013; 5: 669-86.
- Ansermot N, Brawand-Amey M, Kottelat A and Eap CB: Fast quantification of ten psychotropic drugs and metabolites in human plasma by ultra-high performance liquid chromatography tandem mass spectrometry for therapeutic drug monitoring. *J Chromatogr A* 2013; 1292: 160-72.
- Jia M, Li J, He X, Liu M, Zhou Y, Fan Y, *et al.*: Simultaneous determination of iloperidone and its two active metabolites in human plasma by liquid chromatography-tandem mass spectrometry: application to a pharmacokinetic study. *J Chromatogr B Analyt Technol Biomed Life Sci* 2013; 928: 52-7.
- Patteet L, Maudens KE, Stove CP, Lambert WE, Morrens M, Sabbe B, *et al.*: The use of dried blood spots for quantification of 15 antipsychotics and 7 metabolites with ultra-high performance liquid chromatography - tandem mass spectrometry. *Drug Test Anal* 2015; 7: 502-11.
- Mandpe L, Kyadarkunte A and Pokharkar V: High-performance liquid chromatographic method for simultaneous determination of iloperidone and idebenone in spiked plasma. *International Journal of Chemical and Analytical Science* 2013; 4: 19-23.
- Gilby AC: Chromatography: Liquid. Detectors: Ultraviolet and Visible Detection A2 - Wilson, Ian D. *Encyclopedia of Separation Science*. Academic Press, Oxford, 2007: 1-15.
- Es'baghi Z: Photodiode Array Detection in Clinical Applications; Quantitative Analyte Assay Advantages, Limitations and Disadvantages. In: Shi J-W, editor. *Photodiodes - Communications, Bio-Sensings, Measurements and High-Energy Physics*. In Tech, Rijeka, 2011: Ch. 09.
- US Department of Health and Human Services, Food and Drug Administration and Center for Drug Evaluation and Research. *Guidance for Industry, Bioanalytical Method Validation* 2013. Available from: www.fda.gov/downloads/drugs/guidances/ucm368107.pdf.
- Kim S, Thiessen PA, Bolton EE, Chen J, Fu G, Gindulyte A, *et al.*: Pub Chem Substance and Compound databases. *Nucleic Acids Res* 2016; 44: D1202-13.
- Arora K and Gangadharappa HV: An approach to bioanalytical method development and validation: A review. *Int J Pharm Sci Res* 2016; 7: 2291-301.

How to cite this article:

Dubey V and Saini TR: Development and validation of an RP-HPLC method for quantitative estimation of Iloperidone in rat plasma. *Int J Pharm Sci & Res* 2018; 9(2): 790-94. doi: 10.13040/IJPSR.0975-8232.9(2).790-94.